



## Product specification

**V-Monovette® Urine, 9.5 ml, cap yellow, (LxØ): 100 x 15 mm, 50 piece(s)/bag**



## Product information

Product number	11.2552.001
----------------	-------------

## Product descriptions

Name	V-Monovette® Urine, 9.5 ml, Cap yellow, (LxØ): 100 x 15 mm, 50 piece(s)/bag
Product number	11.2552.001
Product description	V-Monovette® Urine, nominal volume: 9.5 ml, without preparation, vacuum system for closed urine transfer, (LxØ): 100 x 15 mm, transparent, fill mark, material: PET, membrane cap, cap yellow, with paper label, sterile, 50 piece(s)/bag

## Material & colours

Product material	Polyethylene therephthalate (PET)
Color of product	transparent
Tube material	Polyethylene therephthalate (PET)
Membrane material	Rubber
Closure material	High Density Polyethylene (HD-PE)
Color of cap	yellow

Diameter	15 mm
Length excluding cap	100 mm
Length including cap	106 mm
Sample volume	9.5 ml

## Product characteristics

Base shape	Conical base
Label/Print	With paper label
Color of print/label	yellow
Graduation	Fill mark
Light protection	no
Cap	Membrane cap
Closure type	Push cap
Type of collection	Vacuum
Centrifugation recommended	400 x g
Centrifugation time	5 Min

## Purity & certification

Product category	diagnostic device, CE
CE certified	CE
Purity standard	sterile
Sterilisation	Gamma irradiation
Batched	yes
Safety data sheet	no

## Packaging

Minimum order quantity	500
Piece(s) / inner box	50
Piece(s) / case	500
Piece(s) / pallet	22000
Height of case	155 mm
Width of case	375 mm
Depth of case	503 mm
Volume of case	0.0292 cbm
Weight of case	4.88 kg
EAN inner box	4038917244886
EAN case	4038917244879

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

This document was prepared by EDP support and is valid without signature.



# SARSTEDT