



Product specification

V-Monovette® Urine, Boric acid, 4 ml, cap green, (LxØ): 75 x 13 mm, 50 piece(s)/bag



Product information

Product number	11.2253.001
----------------	-------------

Product descriptions

Name	V-Monovette® Urine, Boric acid, 4 ml, Cap green, (LxØ): 75 x 13 mm, 50 piece(s)/bag
Product number	11.2253.001
Product description	V-Monovette® Urine, Boric acid, nominal volume: 4 ml, Boric acid, vacuum system for closed urine transfer, (LxØ): 75 x 13 mm, transparent, fill mark, material: PET, membrane cap, cap green, with paper label, sterile, 50 piece(s)/bag

Material & colours

Color of product	transparent
Tube material	Polyethylene therephthalate (PET)
Membrane material	Rubber
Closure material	High Density Polyethylene (HD-PE)
Color of cap	green

Size

Diameter	13 mm
----------	-------

Length excluding cap	75 mm
Length including cap	81 mm
Sample volume	4 ml

Product characteristics

Base shape	Round base
Label/Print	With paper label
Color of print/label	green
Graduation	Fill mark
Light protection	no
Preparation concentration	15 mg/ml
Cap	Membrane cap
Closure type	Push cap
Connection	V-Monovette
Type of collection	Vacuum

Purity & certification

Product category	In-vitro diagnostic (IVD)
CE certified	CE - manufacturer's self-declaration
Purity standard	sterile
Sterilisation	Gamma irradiation
Batched	yes
Safety data sheet	yes

Packaging

Piece(s) / inner box	50
Piece(s) / case	500
Piece(s) / pallet	39000
Height of case	130 mm
Width of case	368 mm
Depth of case	398 mm
Volume of case	0.019 cbm
Weight of case	3.47 kg
EAN inner box	4038917245128
EAN case	4038917245111

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

This document was prepared by EDP support and is valid without signature.



SARSTEDT