




Product specification

Collection tip for Urine Monovette®, 100 piece(s)/bag



Product information

Product number	10.251
----------------	--------

Product descriptions

Name	Collection tip for Urine Monovette®, 100 piece(s)/bag
Product number	10.251
Product description	Collection tip for Urine Monovette®, for urine transfer into Monovette® Urine, 100 piece(s)/bag

Material & colours

Color of product	yellow
------------------	--------

Size

Diameter	8 mm
Length of product	76 mm

Product characteristics

Light protection	no
Connection	Luer (f)
Type of collection	Aspiration

Purity & certification

Product category	In-vitro diagnostic (IVD)
CE certified	CE - manufacturer's self-declaration

Batched	yes
Safety data sheet	no

Packaging

Piece(s) / inner box	100
Piece(s) / case	5000
Piece(s) / pallet	240000
Height of inner box	10 mm
Width of inner box	135 mm
Depth of inner box	200 mm
Height of case	205 mm
Width of case	378 mm
Depth of case	388 mm
Volume of case	0.0301 cbm
Weight of case	3.31 kg
EAN inner box	4038917080941
EAN case	4038917000116

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

This document was prepared by EDP support and is valid without signature.



SARSTEDT