



Product specification

**Monovette® Luer Lithium heparin LH, 4.5 ml, cap green,
(LxØ): 92 x 11 mm, with paper label**



Product information

Product number	05.276.160
----------------	------------

Product descriptions

Name	Monovette® Luer Lithium heparin LH, 4.5 ml, cap green, (LxØ): 92 x 11 mm, with paper label
Product number	05.276.160
Product description	Monovette® Luer Lithium heparin LH, preparation: lithium heparin LH, nominal volume: 4.5 ml, (LxØ): 92 x 11 mm, Luer screw cap, cap green, with paper label, label/print: green, animal venous blood collection, 100 piece(s)/bag

Material & colours

Color of product	transparent
Piston rod material	Polystyrene (PS)
Tube material	Polypropylene (PP)
Piston material	High Density Polyethylene (HD-PE)
Closure material	High Density Polyethylene (HD-PE)
Color of cap	green

Diameter	11 mm
Length excluding cap	92 mm
Length including cap	117 mm
Length of product	92 mm
Height of product	92 mm
Length including cap and plunger	127 mm
Sample volume	4.5 ml

Product characteristics

Label/Print	With paper label
Color of print/label	green
Cap	Luer screw cap
Closure type	Screw cap
Connection	Luer
Type of collection	venous
Centrifugation recommended	2000 x g
Centrifugation time	10 Min

Purity & certification

Product category	no medical device IVD
Batched	yes
Veterinary Use Only	yes

Packaging

Minimum order quantity	100
Piece(s) / inner box	100
Piece(s) / case	500
Piece(s) / pallet	36000
Volume of case	0.0182 cbm
Weight of case	3.6 kg
EAN inner box	4038917104784
EAN case	4038917032162

This document was prepared by EDP support
and is valid without signature.



SARSTEDT