



Product specification

Monovette® Luer Serum CAT, 4.5 ml, cap red, (LxØ): 92 x 11 mm, with paper label



Product information

Product number	05.274.160
----------------	------------

Product descriptions

Name	Monovette® Luer Serum CAT, 4.5 ml, cap red, (LxØ): 92 x 11 mm, with paper label
Product number	05.274.160
Product description	Monovette® Luer Serum CAT, preparation: clotting activator, nominal volume: 4.5 ml, (LxØ): 92 x 11 mm, Luer screw cap, cap red, with paper label, label/print: red, animal venous blood collection, 100 piece(s)/bag

Material & colours

Color of product	transparent
Piston rod material	Polystyrene (PS)
Tube material	Polypropylene (PP)
Piston material	High Density Polyethylene (HD-PE)
Closure material	High Density Polyethylene (HD-PE)
Color of cap	red

Size

Diameter	11 mm
----------	-------

Length excluding cap	92 mm
Length including cap	117 mm
Length of product	92 mm
Height of product	92 mm
Length including cap and plunger	127 mm
Sample volume	4.5 ml

Product characteristics

Label/Print	With paper label
Color of print/label	red
Cap	Luer screw cap
Closure type	Screw cap
Connection	Luer
Type of collection	venous
Centrifugation recommended	2000 x g
Centrifugation time	10 Min

Purity & certification

Product category	no medical device IVD
Batched	yes
Veterinary Use Only	yes

Packaging

Minimum order quantity	100
Piece(s) / inner box	100
Piece(s) / case	500
Piece(s) / pallet	36000
Volume of case	0.0182 cbm
Weight of case	3.6 kg
EAN inner box	4038917104753
EAN case	4038917032322

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

This document was prepared by EDP support
and is valid without signature.



SARSTEDT