




## Product specification

**S-Monovette® Lithium heparin gel LH, 4.7 ml, cap green,  
(LxØ): 75 x 15 mm, with paper label**



## Product information

Product number	03.1631.100
----------------	-------------

## Product descriptions

Name	S-Monovette® Lithium heparin gel LH, 4.7 ml, cap green, (LxØ): 75 x 15 mm, with paper label
Product number	03.1631.100
Product description	S-Monovette® Lithium heparin gel LH, preparation: Lithium heparin gel, 4.7 ml, membrane screw cap, cap green, colour code ISO, (LxØ) without cap: 75 x 15 mm, with paper label, label/print: green, 50 piece(s)/case, sterile

## Material & colours

Color of product	transparent
Piston rod material	Polystyrene (PS)
Tube material	Polypropylene (PP)
Piston material	High Density Polyethylene (HD-PE)
Membrane material	Rubber
Closure material	High Density Polyethylene (HD-PE)
Color of cap	green

## Size

Diameter	15 mm
Length excluding cap	75 mm
Length including cap	91 mm
Length including cap and plunger	112 mm
Sample volume	4.7 ml

## Product characteristics

Label/Print	With paper label
Color of print/label	green
Graduation	Fill mark
Preparation concentration	25 IU/ml blood
Cap	Membrane screw cap
Closure type	Screw cap
Connection	Membrane screw cap
Type of collection	venous
Centrifugation recommended	2500 x g
Centrifugation temperature	20 °C
Centrifugation time	15 Min
Liquid-tight	yes

## Purity & certification

Satisfies the requirement	IATA, ADR
Product category	diagnostic device, CE
CE certified	CE
Purity standard	sterile
Sterilisation	Electron irradiation
Batched	yes

## Packaging

Minimum order quantity	500
Piece(s) / inner box	50
Piece(s) / case	500
Piece(s) / pallet	26000
Height of inner box	82 mm
Width of inner box	117 mm
Depth of inner box	189 mm
Height of case	440 mm
Width of case	129 mm
Depth of case	400 mm

Volume of case	0.0227 cbm
Weight of case	4.54 kg
EAN inner box	4038917132213
EAN case	4038917024327

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

This document was prepared by EDP support and is valid without signature.



# SARSTEDT