



Product specification

**Monovette® Luer Lithium heparin LH, 9 ml, cap natural,
(LxØ): 92 x 16 mm, with paper label**



Product information

Product number	02.265.160
----------------	------------

Product descriptions

Name	Monovette® Luer Lithium heparin LH, 9 ml, cap natural, (LxØ): 92 x 16 mm, with paper label
Product number	02.265.160
Product description	Monovette® Luer Lithium heparin LH, preparation: Lithium heparin, nominal volume: 9 ml, (LxØ): 92 x 16 mm, Luer screw cap, cap natural, with paper label, label/print: green, animal venous blood collection, 100 piece(s)/bag

Material & colours

Color of product	transparent
Piston rod material	Polystyrene (PS)
Tube material	Polypropylene (PP)
Piston material	High Density Polyethylene (HD-PE)
Closure material	High Density Polyethylene (HD-PE)
Color of cap	natural

Size

Diameter	16 mm
----------	-------

Length excluding cap	92 mm
Length including cap	117 mm
Length including cap and plunger	127 mm
Sample volume	9 ml

Product characteristics

Label/Print	With paper label
Color of print/label	green
Cap	Luer screw cap
Closure type	Screw cap
Connection	Luer
Type of collection	venous
Centrifugation recommended	2000 x g
Centrifugation temperature	20 °C
Centrifugation time	10 Min

Purity & certification

Product category	not regulated
CE certified	no CE mark
Batched	yes
Veterinary Use Only	yes

Packaging

Piece(s) / inner box	100
Piece(s) / case	500
Piece(s) / pallet	20000
Volume of case	0.0382 cbm
Weight of case	5.7 kg
EAN inner box	4038917104715
EAN case	4038917032292

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

This document was prepared by EDP support
and is valid without signature.



SARSTEDT