




## Product specification

**S-Monovette® neutral Z, 7.5 ml, cap white, (LxØ): 92 x 15 mm, with paper label**



## Product information

Product number	01.1728.001
----------------	-------------

## Product descriptions

Name	S-Monovette® neutral Z, 7.5 ml, cap white, (LxØ): 92 x 15 mm, with paper label
Product number	01.1728.001
Product description	S-Monovette® neutral Z, without preparation, 7.5 ml, membrane screw cap, cap white, colour code ISO, (LxØ) without cap: 92 x 15 mm, with paper label, label/print: black, 50 piece(s)/case, sterile

## Material & colours

Color of product	transparent
Piston rod material	Polystyrene (PS)
Tube material	Polypropylene (PP)
Piston material	High Density Polyethylene (HD-PE)
Membrane material	Rubber
Closure material	High Density Polyethylene (HD-PE)
Color of cap	white

Diameter	15 mm
Length excluding cap	92 mm
Length including cap	108 mm
Length including cap and plunger	129 mm
Sample volume	7.5 ml

## Product characteristics

Label/Print	With paper label
Color of print/label	black
Graduation	Fill mark
Cap	Membrane screw cap
Closure type	Screw cap
Connection	Membrane screw cap
Type of collection	venous
Centrifugation recommended	2000 x g
Centrifugation temperature	20 °C
Centrifugation time	10 Min
Liquid-tight	yes

## Purity & certification

Satisfies the requirement	IATA, ADR
Product category	In-vitro diagnostic (IVD)
CE certified	CE - manufacturer's self-declaration
Purity standard	sterile
Sterilisation	Electron irradiation
Batched	yes

## Packaging

Piece(s) / inner box	50
Piece(s) / case	500
Piece(s) / pallet	32000
Height of inner box	82 mm
Width of inner box	134 mm
Depth of inner box	189 mm
Height of case	144 mm
Width of case	400 mm
Depth of case	440 mm
Volume of case	0.0253 cbm
Weight of case	4.63 kg
EAN inner box	4038917133326

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

This document was prepared by EDP support and is valid without signature.



# SARSTEDT