




Product specification

PCR strip of 8, PCR Performance Tested, transparent, PP, flat cap



Product information

Product number	72.991.102
----------------	------------

Product descriptions

Name	PCR strip of 8, PCR Performance Tested, transparent, PP, flat cap
Product number	72.991.102
Product description	PCR strip of 8, High Profile, PCR Performance Tested, transparent, material: PP, flat cap, cap attached

Material & colours

Product material	Polypropylene (PP)
Color of product	transparent
Profile	High Profile

Size

Sample volume	0.2 ml
---------------	--------

Product characteristics

Base shape	conical
Graduation	no
Cap	flat cap
Closure type	Snap-on lid

Autoclavable (max. °C)	121 °C
Centrifugation max (RCF)	8000 x g
Centrifugation temperature	20 °C
Centrifugation time	5 Min

Purity & certification

Product category	no medical device IVD
Purity standard	PCR Performance Tested
Batched	yes

Packaging

Minimum order quantity	1200
Piece(s) / inner box	120
Piece(s) / case	1200
Piece(s) / pallet	40320
Volume of case	0.0354 cbm
Weight of case	2.7696 kg
EAN inner box	4038917231657
EAN case	4038917231640

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

This document was prepared by EDP support and is valid without signature.



SARSTEDT