




Product specification

PCR strip of 8, PCR Performance Tested, colour mix, PP, round cap



Product information

| | |
|----------------|------------|
| Product number | 72.994.992 |
|----------------|------------|

Product descriptions

| | |
|---------------------|---------------------------------------------------------------------------------------------------------------------------------------------|
| Name | PCR strip of 8, PCR Performance Tested, Colour mix, PP, round cap |
| Product number | 72.994.992 |
| Product description | PCR strip of 8, High Profile, PCR Performance Tested, red/green/blue/violet/yellow, material: PP, round cap, cap attached, 120 piece(s)/bag |

Material & colours

| | |
|------------------|--------------|
| Color of product | Colour mix |
| Profile | High Profile |

Size

| | |
|---------------|--------|
| Sample volume | 0.2 ml |
|---------------|--------|

Product characteristics

| | |
|------------|-----------|
| Base shape | conical |
| Graduation | no |
| Cap | round cap |

Closure type

Push cap

Purity & certification

| | |
|------------------|------------------------|
| Product category | not regulated |
| CE certified | no CE mark |
| Purity standard | PCR Performance Tested |
| Batched | yes |

Packaging

| | |
|----------------------|---------------|
| Piece(s) / inner box | 120 |
| Piece(s) / case | 1200 |
| Piece(s) / pallet | 1200 |
| Volume of case | 0.0299 cbm |
| Weight of case | 4.2996 kg |
| EAN inner box | 4038917124331 |
| EAN case | 4038917080743 |

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

This document was prepared by EDP support and is valid without signature.



SARSTEDT