




Product specification

PCR strip of 8, PCR Performance Tested, colour mix, PP, round cap



Product information

Product number	72.994.992
----------------	------------

Product descriptions

Name	PCR strip of 8, PCR Performance Tested, Colour mix, PP, round cap
Product number	72.994.992
Product description	PCR strip of 8, High Profile, PCR Performance Tested, red/green/blue/violet/yellow, material: PP, round cap

Material & colours

Product material	Polypropylene (PP)
Color of product	Colour mix
Profile	High Profile

Size

Sample volume	0.2 ml
---------------	--------

Product characteristics

Base shape	conical
Graduation	no
Cap	round cap

Closure type

Push cap

Purity & certification

Product category	no medical device IVD
Purity standard	PCR Performance Tested
Batched	yes

Packaging

Minimum order quantity	1200
Piece(s) / inner box	120
Piece(s) / case	1200
Piece(s) / pallet	1200
Volume of case	0.0299 cbm
Weight of case	4.2996 kg
EAN inner box	4038917124331
EAN case	4038917080743

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

This document was prepared by EDP support and is valid without signature.



SARSTEDT