




# Product specification

Urine container, 90 ml, (ØxH): 60 x 65 mm, PP, with paper label, transparent



## Product information

Product number	75.1435.004
----------------	-------------

## Product descriptions

Name	Urine container, 90 ml, (ØxH): 60 x 65 mm, PP, with paper label, transparent
Product number	75.1435.004
Product description	Urine container, max. working volume: 90 ml, (ØxH): 60 x 65 mm, Ø opening: 65 mm, transparent, graduated, material: PP, screw cap, with paper label, cap assembled, sterile, 400 piece(s)/case

## Material & colours

Product material	Polypropylene (PP)
Color of product	transparent

## Size

Diameter	60 mm
Length excluding cap	60 µm
Height of product	65 mm
Diameter of opening	65 mm
Max. volume	90 ml
Sample volume	90 ml

## Product characteristics

Label/Print	With paper label
Graduation	yes
Closure type	Screw cap

## Purity & certification

Product category	no medical device   IVD
Purity standard	sterile
Sterilisation	Electron irradiation
Batched	yes

## Packaging

Minimum order quantity	400
Piece(s) / inner box	400
Piece(s) / case	400
Piece(s) / pallet	7200
Volume of case	0.1062 cbm
Weight of case	8.23 kg
EAN inner box	4038917385220
EAN case	4038917385213

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

This document was prepared by EDP support and is valid without signature.



# SARSTEDT