



Product specification

Urine container, 75 ml, (ØxH): 65 x 48 mm, PP, transparent



Product information

Product number	75.1437.100
----------------	-------------

Product descriptions

Name	Urine container, 75 ml, (ØxH): 65 x 48 mm, PP, transparent
Product number	75.1437.100
Product description	Urine container, max. working volume: 75 ml, (ØxH): 65 x 48 mm, transparent, graduated, material: PP

Material & colours

Product material	Polypropylene (PP)
Color of product	transparent

Size

Diameter	65 mm
Height of product	48 mm
Max. volume	75 ml
Sample volume	75 ml

Product characteristics

Graduation	yes
------------	-----

Purity & certification

Product category	no medical device IVD
------------------	-------------------------

Batched

yes

Packaging

Minimum order quantity	1000
Piece(s) / inner box	25
Piece(s) / case	1000
Piece(s) / pallet	28000
Volume of case	0.065 cbm
Weight of case	9.65 kg
EAN inner box	4038917129138
EAN case	4038917074711

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

This document was prepared by EDP support and is valid without signature.



SARSTEDT