



## Product specification

**Urine container, 148 ml, (ØxH): 73 x 67 mm, PS, white**



### Product information

Product number	75.560.901
----------------	------------

### Product descriptions

Name	Urine container, 148 ml, (ØxH): 73 x 67 mm, PS, white
Product number	75.560.901
Product description	Urine container, max. working volume: 148 ml, (ØxH): 73 x 67 mm, white, material: PS

### Material & colours

Product material	Polystyrene (PS)
Color of product	white

### Size

Diameter	73 mm
Height of product	67 mm
Sample volume	148 ml

### Product characteristics

Graduation	no
------------	----

### Purity & certification

Product category	no medical device   IVD
Batched	yes

## Packaging

Minimum order quantity	2500
Piece(s) / inner box	100
Piece(s) / case	2500
Piece(s) / pallet	90000
Height of case	356 mm
Width of case	368 mm
Depth of case	457 mm
Volume of case	0.0599 cbm
Weight of case	8.16 kg
EAN inner box	4038917318198
EAN case	4038917916547

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

This document was prepared by EDP support and is valid without signature.



# SARSTEDT