




Product specification

Rotor, for 36 collection tubes, for micro tube Ø 10.8 mm



Product information

Product number	92.181.610
----------------	------------

Product descriptions

Name	Rotor, for 36 Collection tubes, for micro tube Ø 10.8 mm
Product number	92.181.610
Product description	Rotor, for 36 collection tubes, for micro tube Ø 10.8 mm

Product characteristics

Storage capacity	36
------------------	----

Purity & certification

Product category	not regulated
CE certified	no CE mark
Batched	yes

Packaging

Piece(s) / pallet	40
Volume of case	0.0056 cbm
Weight of case	1.36 kg
EAN inner box	4038917156233
EAN case	4038917156226

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

This document was prepared by EDP support and is valid without signature.



SARSTEDT