



Product specification

Rotor, for 22 collection tubes



Product information

Product number	92.181.611
----------------	------------

Product descriptions

Name	Rotor, for 22 Collection tubes
Product number	92.181.611
Product description	Rotor, for 22 collection tubes

Product characteristics

Storage capacity	22
------------------	----

Purity & certification

Product category	not regulated
CE certified	no CE mark
Batched	yes

Packaging

Piece(s) / pallet	40
Volume of case	0.0052 cbm
Weight of case	1.36 kg
EAN inner box	4038917156257
EAN case	4038917156240

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

This document was prepared by EDP support and is valid without signature.



SARSTEDT