




## Product specification

**Film, free from DNase/RNase, transparent, pierceable**



### Product information

Product number	95.2092
----------------	---------

### Product descriptions

Name	Film, free from DNase/RNase, transparent, pierceable
Product number	95.2092
Product description	Film, free from DNase/RNase, transparent, suitable for qPCR, pierceable, application area: -80 - 120 °C

### Material & colours

Color of product	transparent
------------------	-------------

### Size

Width of product	82 mm
Length of product	125 mm

### Product characteristics

Light protection	no
Closure type	Film
Pierceable	yes
Suitable for qPCR	yes
Minimum temperature	-80 °C
Maximum temperature	120 °C

## Purity & certification

Product category	not regulated
CE certified	no CE mark
Purity standard	free from DNase/RNase
Batched	yes

## Packaging

Piece(s) / inner box	100
Piece(s) / case	1200
Piece(s) / pallet	57600
Height of inner box	35 mm
Width of inner box	135 mm
Depth of inner box	135 mm
Volume of case	0.0142 cbm
Weight of case	4.3596 kg
EAN inner box	4038917335676
EAN case	4038917335669

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

This document was prepared by EDP support and is valid without signature.



# SARSTEDT