




Product specification

Film, free from DNase/RNase, transparent, pierceable



Product information

Product number	95.2093
----------------	---------

Product descriptions

Name	Film, free from DNase/RNase, transparent, pierceable
Product number	95.2093
Product description	Film, free from DNase/RNase, transparent, pierceable, application area: -80 - 120 °C

Material & colours

Color of product	transparent
------------------	-------------

Size

Width of product	75 mm
Length of product	118 mm

Product characteristics

Light protection	no
Closure type	Film
Pierceable	yes
Suitable for qPCR	no
Minimum temperature	-80 °C
Maximum temperature	120 °C

Purity & certification

Product category	not regulated
CE certified	no CE mark
Purity standard	free from DNase/RNase
Batched	yes

Packaging

Piece(s) / inner box	20
Piece(s) / case	240
Piece(s) / pallet	240
Height of inner box	27 mm
Width of inner box	135 mm
Depth of inner box	135 mm
Volume of case	0.0142 cbm
Weight of case	3 kg
EAN inner box	4038917335706
EAN case	4038917335690

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

This document was prepared by EDP support and is valid without signature.



SARSTEDT