



Product specification

**S-Monovette® Lithium heparin LH, 7.5 ml, cap green,
(LxØ): 92 x 15 mm, with paper label**



Product information

Product number	01.1604.100
----------------	-------------

Product descriptions

Name	S-Monovette® Lithium heparin LH, 7.5 ml, cap green, (LxØ): 92 x 15 mm, with paper label
------	-----------------------------------------------------------------------------------------

Product number	01.1604.100
----------------	-------------

Product description	S-Monovette® Lithium heparin LH, preparation: lithium heparin LH, 7.5 ml, membrane screw cap, cap green, colour code ISO, (LxØ) without cap: 92 x 15 mm, with paper label, label/print: green, 50 piece(s)/case, sterile
---------------------	--------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------

Material & colours

Color of product	transparent
Piston rod material	Polystyrene (PS)
Tube material	Polypropylene (PP)
Piston material	High Density Polyethylene (HD-PE)
Membrane material	Rubber
Closure material	High Density Polyethylene (HD-PE)
Color of cap	green

Size

Diameter	15 mm
Length excluding cap	92 mm
Length including cap	108 mm
Length including cap and plunger	129 mm
Sample volume	7.5 ml

Product characteristics

Label/Print	With paper label
Color of print/label	green
Graduation	Fill mark
Preparation concentration	16 IU/ml blood
Cap	Membrane screw cap
Closure type	Screw cap
Connection	Membrane screw cap
Type of collection	venous
Centrifugation recommended	2000 x g
Centrifugation temperature	20 °C
Centrifugation time	10 Min
Liquid-tight	yes

Purity & certification

Satisfies the requirement	IATA, ADR
Product category	diagnostic device, CE
CE certified	CE
Purity standard	sterile
Sterilisation	Electron irradiation
Batched	yes

Packaging

Minimum order quantity	500
Piece(s) / inner box	50
Piece(s) / case	500
Piece(s) / pallet	32000
Height of inner box	82 mm
Width of inner box	134 mm
Depth of inner box	189 mm
Height of case	440 mm
Width of case	144 mm
Depth of case	400 mm

Volume of case	0.0253 cbm
Weight of case	5.05 kg
EAN inner box	4038917131919
EAN case	4038917024037

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

This document was prepared by EDP support and is valid without signature.



SARSTEDT