



Product specification

S-Monovette® EDTA K2E, 2.7 ml, cap violet, (LxØ): 75 x 13 mm, with paper label



Product information

Product number	04.1915.100
----------------	-------------

Product descriptions

Name	S-Monovette® EDTA K2E, 2.7 ml, cap violet, (LxØ): 75 x 13 mm, with paper label
Product number	04.1915.100
Product description	S-Monovette® EDTA K2E, preparation: K2 EDTA, 2.7 ml, membrane screw cap, cap violet, colour code ISO, (LxØ) without cap: 75 x 13 mm, with paper label, label/print: violet, 50 piece(s)/case, sterile

Material & colours

Color of product	transparent
Piston rod material	Polystyrene (PS)
Tube material	Polypropylene (PP)
Piston material	High Density Polyethylene (HD-PE)
Membrane material	Rubber
Closure material	High Density Polyethylene (HD-PE)
Color of cap	violet

Size

Diameter	13 mm
Length excluding cap	75 mm
Length including cap	91 mm
Length including cap and plunger	112 mm
Sample volume	2.7 ml

Product characteristics

Label/Print	With paper label
Color of print/label	violet
Graduation	Fill mark
Preparation concentration	1.6 mg/ml blood
Cap	Membrane screw cap
Closure type	Screw cap
Connection	Membrane screw cap
Type of collection	venous
Liquid-tight	yes

Purity & certification

Satisfies the requirement	IATA, ADR
Product category	diagnostic device, CE
CE certified	CE
Purity standard	sterile
Sterilisation	Electron irradiation
Batched	yes

Packaging

Minimum order quantity	500
Piece(s) / inner box	50
Piece(s) / case	500
Piece(s) / pallet	39000
Height of inner box	73 mm
Width of inner box	116 mm
Depth of inner box	171 mm
Height of case	135 mm
Width of case	358 mm
Depth of case	388 mm
Volume of case	0.0188 cbm
Weight of case	2.86 kg
EAN inner box	4038917382250

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

This document was prepared by EDP support and is valid without signature.



SARSTEDT