



Product specification

**Blood Gas Monovette®, calcium-balanced lithium heparin,
2 ml, cap white/orange, connection: Luer (m)**



Product information

Product number	05.1147.028
----------------	-------------

Product descriptions

Name	Blood Gas Monovette®, calcium-balanced lithium heparin, 2 ml, cap white/orange, connection: Luer (m)
------	--

Product number	05.1147.028
----------------	-------------

Product description	Blood Gas Monovette®, preparation: calcium-balanced lithium heparin, nominal volume: 2 ml, cap white/orange, label/print: transparent/orange, connection: Luer (m), (LxØ) without cap: 66 x 11 mm, with membrane adapter assembled, 1 piece(s)/blister, sterile, 100 piece(s)/case
---------------------	--

Material & colours

Color of product	transparent
------------------	-------------

Piston rod material	Polystyrene (PS)
---------------------	------------------

Tube material	Polypropylene (PP)
---------------	--------------------

Piston material	High Density Polyethylene (HD-PE)
-----------------	-----------------------------------

Closure material	High Density Polyethylene (HD-PE)
------------------	-----------------------------------

Color of cap	white/orange
--------------	--------------

Size

Diameter	11 mm
Length excluding cap	66 mm
Length including cap	82 mm
Length including cap and plunger	100 mm
Sample volume	2 ml

Product characteristics

Design	Membrane adapter assembled
Label/Print	With plastic label
Color of print/label	transparent/orange
Graduation	Fill mark
Preparation concentration	50 IU/ml blood
Connection	Luer (m)
Type of collection	Arterial and venous

Purity & certification

Product category	In-vitro diagnostic (IVD)
CE certified	CE - manufacturer's self-declaration
Purity standard	sterile
Sterilisation	Electron irradiation
Batched	yes

Packaging

Type of inner packaging	Case
Piece(s) / inner box	100
Piece(s) / case	500
Piece(s) / pallet	18000
Height of inner box	171 mm
Width of inner box	111 mm
Depth of inner box	351 mm
Height of case	195 mm
Width of case	368 mm
Depth of case	578 mm
Volume of case	0.0415 cbm
Weight of case	3.4 kg
EAN inner box	4038917280679
EAN case	4038917280662

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

This document was prepared by EDP support and is valid without signature.



SARSTEDT