

Product specification

S-Monovette® K2 EDTA, 1.2 ml, Cap violet, (LxØ): 66 x 8 mm, with paper label



Product information

Order number	06.1660.100
--------------	-------------

Product descriptions

Name	S-Monovette® K2 EDTA, 1.2 ml, Cap violet, (LxØ): 66 x 8 mm, with paper label
Product number	06.1660.100
Product description	S-Monovette® 1.2 ml, K2 EDTA, 66x8 mm, violet ISO code, transparent label, 50/inner box sterile

Product characteristics

Label/Print	With paper label
Color of print/label	violet
Graduation	yes
Preparation concentration	1.6 mg/ml blood
Cap	Membrane screw cap
Type of cap	Screw cap
Connection	Membrane screw cap
Type of collection	venous
Liquid-tight	yes

Size

Diameter	8 mm
Length excluding cap	66 mm
Length including cap	82 mm
Length including cap and plunger	103 mm
Sample volume	1.2 ml

Material & colours

Color of product	transparent
Piston rod material	Polystyrene (PS)
Tube material	Polypropylene (PP)
Piston material	High Density Polyethylene (HD-PE)
Membrane material	Rubber
Cap material	High Density Polyethylene (HD-PE)
Color of cap	violet

Purity & certification

Satisfies the requirement	IATA, ADR
Accordance with MDD or IVDD	IVD product
CE certified	CE
Purity standard	sterile
Sterilisation	Electron irradiation
Lot or serial number	yes

Packaging

Piece(s) / inner box	50
Piece(s) / case	500
Piece(s) / pallet	56000
Height of inner box	63 mm
Width of inner box	108 mm
Depth of inner box	145 mm
Height of case	125 mm
Width of case	303 mm
Depth of case	333 mm
Volume of case	0.0126 cbm
Weight of case	2.18 kg
EAN inner box	4038917322737
EAN case	4038917322720

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

This document was prepared by EDP support and is valid without signature.

