



## Product specification

**S-Monovette® Serum Gel CAT, 1.1 ml, cap brown, (LxØ):  
66 x 8 mm, with plastic label**



## Product information

|                |             |
|----------------|-------------|
| Product number | 06.1667.001 |
|----------------|-------------|

## Product descriptions

|                     |   |
|---------------------|---|
| Name                | S-Monovette® Serum Gel CAT, 1.1 ml, cap brown, (LxØ): 66 x 8 mm, with plastic label   |
| Product number      | 06.1667.001   |
| Product description | S-Monovette® Serum Gel CAT, preparation: clotting activator/gel, 1.1 ml, membrane screw cap, cap brown, colour code EU/ISO, (LxØ) without cap: 66 x 8 mm, with plastic label, label/print: white/brown, 50 piece(s)/case, sterile |

## Material & colours

|                     |                                   |
|---------------------|-----------------------------------|
| Color of product    | transparent                       |
| Piston rod material | Polystyrene (PS)                  |
| Tube material       | Polypropylene (PP)                |
| Piston material     | rubber                            |
| Membrane material   | Rubber                            |
| Closure material    | High Density Polyethylene (HD-PE) |
| Color of cap        | brown                             |

## Size

|                                  |        |
|----------------------------------|--------|
| Diameter                         | 8 mm   |
| Length excluding cap             | 66 mm  |
| Length including cap             | 82 mm  |
| Length including cap and plunger | 103 mm |
| Sample volume                    | 1.1 ml |

## Product characteristics

|                            |                    |
|----------------------------|--------------------|
| Label/Print                | With plastic label |
| Color of print/label       | white/brown        |
| Graduation                 | yes                |
| Cap                        | Membrane screw cap |
| Closure type               | Screw cap          |
| Connection                 | Membrane screw cap |
| Type of collection         | venous             |
| Centrifugation recommended | 2500 x g           |
| Centrifugation temperature | 20 °C              |
| Centrifugation time        | 10 Min             |
| Liquid-tight               | yes                |

## Purity & certification

|                           |                                      |
|---------------------------|--------------------------------------|
| Satisfies the requirement | IATA, ADR                            |
| Product category          | In-vitro diagnostic (IVD)            |
| CE certified              | CE - manufacturer's self-declaration |
| Purity standard           | sterile                              |
| Sterilisation             | Electron irradiation                 |
| Batched                   | yes                                  |

## Packaging

|                      |            |
|----------------------|------------|
| Piece(s) / inner box | 50         |
| Piece(s) / case      | 500        |
| Piece(s) / pallet    | 56000      |
| Height of inner box  | 63 mm      |
| Width of inner box   | 108 mm     |
| Depth of inner box   | 145 mm     |
| Height of case       | 125 mm     |
| Width of case        | 303 mm     |
| Depth of case        | 333 mm     |
| Volume of case       | 0.0126 cbm |
| Weight of case       | 2.66 kg    |

EAN inner box

4038917132893

EAN case

4038917024532

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

This document was prepared by EDP support and is valid without signature.



# SARSTEDT