



## Product specification

**Urine Monovette®, 10 ml, cap yellow, (LxØ): 102 x 15 mm,  
1 piece(s)/blister**



## Product information

Product number	10.252.020
----------------	------------

## Product descriptions

Name	Urine Monovette®, 10 ml, cap yellow, (LxØ): 102 x 15 mm, 1 piece(s)/blister
------	--------------------------------------------------------------------------------

Product number	10.252.020
----------------	------------

Product description	Urine Monovette®, nominal volume: 10 ml, without preparation, Urine tube, (LxØ): 102 x 15 mm, transparent, graduated, material: PP, Luer screw cap, cap yellow, incl. collection tip, with plastic label, sterile, 1 piece(s)/blister
---------------------	---------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------

## Material & colours

Color of product	transparent
------------------	-------------

Piston rod material	Polystyrene (PS)
---------------------	------------------

Tube material	Polypropylene (PP)
---------------	--------------------

Piston material	High Density Polyethylene (HD-PE)
-----------------	-----------------------------------

Closure material	High Density Polyethylene (HD-PE)
------------------	-----------------------------------

Color of cap	yellow
--------------	--------

## Size

Diameter	15 mm
----------	-------

Length excluding cap	102 mm
Length including cap	117 mm
Length including cap and plunger	127 mm
Sample volume	10 ml

## Product characteristics

Base shape	Round base
Label/Print	With plastic label
Color of print/label	transparent/yellow
Graduation	yes
Light protection	no
Cap	Luer screw cap
Closure type	Screw cap
Connection	Luer (m)
Example of use	Urine collection, transport, and storage
Type of collection	Aspiration
Centrifugation recommended	400 x g
Centrifugation time	5 Min

## Purity & certification

Satisfies the requirement	IATA, ADR
Product category	In-vitro diagnostic (IVD)
CE certified	CE - manufacturer's self-declaration
Purity standard	sterile
Sterilisation	Electron irradiation
Batched	yes

## Packaging

Piece(s) / inner box	100
Piece(s) / case	500
Piece(s) / pallet	13500
Height of inner box	171 mm
Width of inner box	151 mm
Depth of inner box	351 mm
Height of case	195 mm
Width of case	368 mm
Depth of case	778 mm
Volume of case	0.0558 cbm
Weight of case	5.82 kg
EAN inner box	4038917104432

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

This document was prepared by EDP support and is valid without signature.



# SARSTEDT