



## Product specification

**V-Monovette® Urine, 10 ml, cap yellow, (LxØ): 100 x 15 mm, 50 piece(s)/bag**



## Product information

|                |             |
|----------------|-------------|
| Product number | 11.2452.001 |
|----------------|-------------|

## Product descriptions

|                     |  |
|---------------------|--|
| Name                | V-Monovette® Urine, 10 ml, Cap yellow, (LxØ): 100 x 15 mm, 50 piece(s)/bag   |
| Product number      | 11.2452.001  |
| Product description | V-Monovette® Urine, nominal volume: 10 ml, without preparation, vacuum system for closed urine transfer, (LxØ): 100 x 15 mm, transparent, fill mark, material: PET, membrane cap, cap yellow, with paper label, sterile, 50 piece(s)/bag |

## Material & colours

|                   |                                   |
|-------------------|-----------------------------------|
| Color of product  | transparent                       |
| Tube material     | Polyethylene therephthalate (PET) |
| Membrane material | Rubber                            |
| Closure material  | High Density Polyethylene (HD-PE) |
| Color of cap      | yellow                            |

## Size

|          |       |
|----------|-------|
| Diameter | 15 mm |
|----------|-------|

|                      |        |
|----------------------|--------|
| Length excluding cap | 100 mm |
| Length including cap | 106 mm |
| Sample volume        | 10 ml  |

## Product characteristics

|                            |                  |
|----------------------------|------------------|
| Base shape                 | Round base       |
| Label/Print                | With paper label |
| Color of print/label       | yellow           |
| Graduation                 | Fill mark        |
| Light protection           | no               |
| Cap                        | Membrane cap     |
| Closure type               | Push cap         |
| Type of collection         | Vacuum           |
| Centrifugation recommended | 400 x g          |
| Centrifugation time        | 5 Min            |

## Purity & certification

|                   |                                      |
|-------------------|--------------------------------------|
| Product category  | In-vitro diagnostic (IVD)            |
| CE certified      | CE - manufacturer's self-declaration |
| Purity standard   | sterile                              |
| Sterilisation     | Gamma irradiation                    |
| Batched           | yes                                  |
| Safety data sheet | no                                   |

## Packaging

|                      |               |
|----------------------|---------------|
| Piece(s) / inner box | 50            |
| Piece(s) / case      | 500           |
| Piece(s) / pallet    | 22000         |
| Height of case       | 155 mm        |
| Width of case        | 375 mm        |
| Depth of case        | 503 mm        |
| Volume of case       | 0.0292 cbm    |
| Weight of case       | 4.95 kg       |
| EAN inner box        | 4038917244855 |
| EAN case             | 4038917531948 |

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

This document was prepared by EDP support and is valid without signature.



# SARSTEDT