



Product specification

Multi-adapter, With pin



Product information

Product number	14.1205.300
----------------	-------------

Product descriptions

Name	Multi-adapter, With pin
Product number	14.1205.300
Product description	Multi-adapter, With pin, connection: S-Monovette®, 1 piece(s)/blister, 100 piece(s)/case, sterile, non-pyrogenic/endotoxin-free

Material & colours

Product material	Polyethylene therephthalate (PET)
Color of product	translucent

Size

Length of product	33 mm
-------------------	-------

Product characteristics

Design	With pin
Connection	S-Monovette

Purity & certification

Product category	no medical device IVD
Purity standard	sterile, non-pyrogenic/endotoxin-free
Sterilisation	Gamma irradiation

Pyrogen-free	yes
Batched	yes

Packaging

Minimum order quantity	50000
Piece(s) / inner box	100
Piece(s) / case	1000
Piece(s) / pallet	80000
Height of inner box	71 mm
Width of inner box	81 mm
Depth of inner box	268 mm
Height of case	287 mm
Width of case	178 mm
Depth of case	373 mm
Volume of case	0.0191 cbm
Weight of case	4.94 kg
EAN inner box	4038917303958
EAN case	4038917303941

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

This document was prepared by EDP support and is valid without signature.



SARSTEDT