



Product specification

Colour-coded cap, pink



Product information

Product number	14.1513
----------------	---------

Product descriptions

Name	Colour-coded cap, pink
Product number	14.1513
Product description	Colour-coded cap, for S-Monovette®, pink, 100 piece(s)/bag

Material & colours

Product material	High Density Polyethylene (HD-PE)
Color of product	pink

Size

Length of product	13 mm
-------------------	-------

Product characteristics

Connection	S-Monovette
------------	-------------

Purity & certification

Product category	no medical device   IVD
Batched	no

Packaging

Minimum order quantity	1000
------------------------	------

Piece(s) / inner box	100
Piece(s) / case	1000
Piece(s) / pallet	324000
Height of case	100 mm
Width of case	191 mm
Depth of case	268 mm
Volume of case	0.0051 cbm
Weight of case	0.47 kg
EAN inner box	4038917086790
EAN case	4038917073943

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

This document was prepared by EDP support and is valid without signature.



# SARSTEDT