



## Product specification

### Colour-coded cap, yellow



## Product information

Product number	14.1513.002
----------------	-------------

## Product descriptions

Name	Colour-coded cap, yellow
Product number	14.1513.002
Product description	Colour-coded cap, for S-Monovette®, yellow, 100 piece(s)/bag

## Material & colours

Color of product	yellow
------------------	--------

## Size

Length of product	13 mm
-------------------	-------

## Purity & certification

Product category	not regulated
CE certified	no CE mark
Batched	no

## Packaging

Piece(s) / inner box	100
Piece(s) / case	1000
Piece(s) / pallet	324000
Height of case	100 mm

Width of case	191 mm
Depth of case	268 mm
Volume of case	0.0051 cbm
Weight of case	0.48 kg
EAN inner box	4038917251556
EAN case	4038917251549

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

This document was prepared by EDP support and is valid without signature.



# SARSTEDT