



## Product specification

### Multivette® 600 Serum CAT, 600 µl, cap red, screw cap



## Product information

Product number	15.1670.100
----------------	-------------

## Product descriptions

Name	Multivette® 600 Serum CAT, 600 µl, cap red, screw cap
Product number	15.1670.100
Product description	Multivette® 600 Serum CAT, venous and capillary blood samples, preparation: clotting activator, nominal volume: 600 µl, (LxØ) with cap: 47.9 x 8 mm, with print, cap red, screw cap, flat base, 100 piece(s)/case

## Material & colours

Product material	Polypropylene (PP)
Color of product	transparent
Closure material	High Density Polyethylene (HD-PE)
Color of cap	red

## Size

Diameter	8 mm
Length excluding cap	44.4 mm
Length including cap	47.9 mm
Sample volume	600 µl

## Product characteristics

Base shape	Flat base
Label/Print	With print
Color of print/label	red
Closure type	Screw cap
Type of collection	Venous and capillary
Centrifugation recommended	10000 U/min
Centrifugation time	5 Min

## Purity & certification

Product category	diagnostic device, CE
CE certified	CE
Batched	yes

## Packaging

Minimum order quantity	1000
Piece(s) / inner box	100
Piece(s) / case	1000
Piece(s) / pallet	63000
Height of inner box	81 mm
Width of inner box	81 mm
Depth of inner box	201 mm
Height of case	185 mm
Width of case	213 mm
Depth of case	423 mm
Volume of case	0.0167 cbm
Weight of case	2.79 kg
EAN inner box	4038917132947
EAN case	4038917024587

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.



# SARSTEDT

This document was prepared by EDP support and is valid without signature.