



Product specification

Multivette® 600 EDTA K3E, 600 µl, cap violet, screw cap



Product information

Product number	15.1671.100
----------------	-------------

Product descriptions

Name	Multivette® 600 EDTA K3E, 600 µl, cap violet, screw cap
Product number	15.1671.100
Product description	Multivette® 600 EDTA K3E, venous and capillary blood samples, preparation: K3 EDTA, nominal volume: 600 µl, (LxØ) with cap: 47.9 x 8 mm, with print, cap violet, screw cap, flat base, 100 piece(s)/case

Material & colours

Color of product	transparent
Closure material	High Density Polyethylene (HD-PE)
Color of cap	violet

Size

Diameter	8 mm
Length excluding cap	44.4 mm
Length including cap	47.9 mm
Sample volume	600 µl

Product characteristics

Base shape	Flat base
------------	-----------

Label/Print	With print
Color of print/label	violet
Closure type	Screw cap
Type of collection	Venous and capillary

Purity & certification

Product category	In-vitro diagnostic (IVD)
CE certified	CE - manufacturer's self-declaration
Batched	yes

Packaging

Piece(s) / inner box	100
Piece(s) / case	1000
Piece(s) / pallet	63000
Height of inner box	81 mm
Width of inner box	81 mm
Depth of inner box	201 mm
Height of case	185 mm
Width of case	213 mm
Depth of case	423 mm
Volume of case	0.0167 cbm
Weight of case	2.79 kg
EAN inner box	4038917132954
EAN case	4038917024594

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

This document was prepared by EDP support and is valid without signature.



SARSTEDT