



Product specification

Multivette® 600 Serum Gel CAT, 600 µl, cap brown, screw cap



Product information

Product number	15.1674
----------------	---------

Product descriptions

Name	Multivette® 600 Serum Gel CAT, 600 µl, cap brown, screw cap
------	---

Product number	15.1674
----------------	---------

Product description	Multivette® 600 Serum Gel CAT, venous and capillary blood samples, preparation: clotting activator/gel, nominal volume: 600 µl, (LxØ) with cap: 53.3 x 8 mm, with print, cap brown, screw cap, flat base, 100 piece(s)/case
---------------------	---

Material & colours

Product material	Polypropylene (PP)
Color of product	transparent
Closure material	High Density Polyethylene (HD-PE)
Color of cap	brown

Size

Diameter	8 mm
Length excluding cap	49.8 mm
Length including cap	53.3 mm

Sample volume	600 µl
---------------	--------

Product characteristics

Base shape	Flat base
Label/Print	With print
Color of print/label	brown
Cap	standard screw cap
Closure type	Screw cap
Type of collection	Venous and capillary
Centrifugation recommended	10000 x g
Centrifugation time	5 Min

Purity & certification

Product category	diagnostic device, CE
CE certified	CE
Batched	yes

Packaging

Minimum order quantity	1000
Piece(s) / inner box	100
Piece(s) / case	1000
Piece(s) / pallet	63000
Height of inner box	81 mm
Width of inner box	81 mm
Depth of inner box	201 mm
Height of case	185 mm
Width of case	213 mm
Depth of case	423 mm
Volume of case	0.0167 cbm
Weight of case	2.97 kg
EAN inner box	4038917087421
EAN case	4038917005852

This document was prepared by EDP support
and is valid without signature.



SARSTEDT