



## Product specification

**Minivette® POCT neutral Z, 100 µl, plunger white, 200 piece(s)/bag**



### Product information

Product number	17.2111.100
----------------	-------------

### Product descriptions

Name	Minivette® POCT neutral Z, 100 µl, plunger white, 200 piece(s)/bag
Product number	17.2111.100
Product description	Minivette® POCT neutral Z, nominal volume: 100 µl, without preparation, white, 200 piece(s)/bag, 200 piece(s)/case

### Material & colours

Color of product	transparent
------------------	-------------

### Size

Sample volume	100 µl
---------------	--------

### Product characteristics

Type of collection	capillary
--------------------	-----------

### Purity & certification

Product category	In-vitro diagnostic (IVD)
CE certified	CE - manufacturer's self-declaration

Latex-free	yes
Batched	yes

## Packaging

Piece(s) / inner box	200
Piece(s) / case	2000
Piece(s) / pallet	144000
Height of inner box	117 mm
Width of inner box	99 mm
Depth of inner box	140 mm
Height of case	140 mm
Width of case	293 mm
Depth of case	513 mm
Volume of case	0.021 cbm
Weight of case	2.47 kg
EAN inner box	4038917339278
EAN case	4038917339261

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

This document was prepared by EDP support and is valid without signature.



# SARSTEDT