



Product specification

Minivette® POCT neutral Z, 200 µl, plunger white, 150 piece(s)/bag



Product information

Product number	17.2111.200
----------------	-------------

Product descriptions

Name	Minivette® POCT neutral Z, 200 µl, plunger white, 150 piece(s)/bag
Product number	17.2111.200
Product description	Minivette® POCT neutral Z, nominal volume: 200 µl, without preparation, white, 150 piece(s)/bag, 150 piece(s)/case

Material & colours

Product material	Cellulose acetate (CA)
Color of product	transparent

Size

Sample volume	200 µl
---------------	--------

Product characteristics

Type of collection	capillary
--------------------	-----------

Purity & certification

Product category	diagnostic device, CE
------------------	-----------------------

CE certified	CE
Latex-free	yes
Batched	yes

Packaging

Minimum order quantity	1500
Piece(s) / inner box	150
Piece(s) / case	1500
Piece(s) / pallet	108000
Height of inner box	117 mm
Width of inner box	99 mm
Depth of inner box	140 mm
Height of case	140 mm
Width of case	293 mm
Depth of case	513 mm
Volume of case	0.021 cbm
Weight of case	2.34 kg
EAN inner box	4038917344005
EAN case	4038917343992

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

This document was prepared by EDP support and is valid without signature.



SARSTEDT