



Product specification

Minivette® POCT Lithium heparin LH, 50 µl, plunger green, colour code ISO, 200 piece(s)/bag



Product information

Product number	17.2112.150
----------------	-------------

Product descriptions

Name	Minivette® POCT Lithium heparin LH, 50 µl, plunger green, colour code ISO, 200 piece(s)/bag
------	---

Product number	17.2112.150
----------------	-------------

Product description	Minivette® POCT Lithium heparin LH, nominal volume: 50 µl, preparation: Heparin, green, colour code ISO, 200 piece(s)/bag, 200 piece(s)/case
---------------------	--

Material & colours

Color of product	transparent
------------------	-------------

Size

Sample volume	50 µl
---------------	-------

Product characteristics

Type of collection	capillary
--------------------	-----------

Purity & certification

Product category	In-vitro diagnostic (IVD)
------------------	---------------------------

CE certified	CE - manufacturer's self-declaration
--------------	--------------------------------------

Latex-free	yes
Batched	yes

Packaging

Piece(s) / inner box	200
Piece(s) / case	2000
Piece(s) / pallet	144000
Height of inner box	117 mm
Width of inner box	99 mm
Depth of inner box	140 mm
Height of case	140 mm
Width of case	293 mm
Depth of case	513 mm
Volume of case	0.021 cbm
Weight of case	2.68 kg
EAN inner box	4038917320023
EAN case	4038917320016

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

This document was prepared by EDP support and is valid without signature.



SARSTEDT