



## Product specification

### Microvette® CB 200 EDTA K3E, 200 µl, cap violet, push cap, flat base



## Product information

Product number	18.1321.102
----------------	-------------

## Product descriptions

Name	Microvette® CB 200 EDTA K3E, 200 µl, cap violet, push cap, flat base
------	--

Product number	18.1321.102
----------------	-------------

Product description	Microvette® CB 200 EDTA K3E, capillary blood collection, preparation: K3 EDTA, nominal volume: 200 µl, (LxØ) with cap: 46 x 10.8 mm, with plastic label, cap violet, push cap, flat base, 50 piece(s)/bag
---------------------	---

## Material & colours

Color of product	transparent
Closure material	Low Density Polyethylene (LD-PE)
Color of cap	violet

## Size

Diameter	10.8 mm
Length excluding cap	40 mm
Length including cap	46 mm
Length of product	40 mm
Sample volume	200 µl

## Product characteristics

Base shape	Flat base
Label/Print	With plastic label
Closure type	Push cap
Type of collection	capillary

## Purity & certification

Product category	In-vitro diagnostic (IVD)
CE certified	CE - manufacturer's self-declaration
Batched	yes

## Packaging

Piece(s) / inner box	50
Piece(s) / case	1000
Piece(s) / pallet	48000
Height of case	205 mm
Width of case	378 mm
Depth of case	388 mm
Volume of case	0.0301 cbm
Weight of case	2.51 kg
EAN inner box	4038917345842
EAN case	4038917345835

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

This document was prepared by EDP support and is valid without signature.



# SARSTEDT