




Product specification

Microvette® CB 200 Serum Gel CAT, 200 µl, cap brown, push cap, flat base



Product information

Product number	18.1324
----------------	---------

Product descriptions

Name	Microvette® CB 200 Serum Gel CAT, 200 µl, cap brown, push cap, flat base
Product number	18.1324
Product description	Microvette® CB 200 Serum Gel CAT, capillary blood collection, nominal volume: 200 µl, (LxØ) with cap: 46 x 10.8 mm, cap brown, push cap, flat base, 100 piece(s)/bag

Material & colours

Product material	Polypropylene (PP)
Color of product	transparent
Closure material	Low Density Polyethylene (LD-PE)
Color of cap	brown

Size

Diameter	10.8 mm
Length excluding cap	40 mm
Length including cap	46 mm
Length of product	40 mm

Sample volume	200 µl
---------------	--------

Product characteristics

Base shape	Flat base
Closure type	Push cap
Type of collection	capillary

Purity & certification

Product category	diagnostic device, CE
CE certified	CE
Batched	yes

Packaging

Minimum order quantity	1000
Piece(s) / inner box	100
Piece(s) / case	1000
Piece(s) / pallet	64000
Height of case	215 mm
Width of case	298 mm
Depth of case	398 mm
Volume of case	0.0255 cbm
Weight of case	2.46 kg
EAN inner box	4038917161770
EAN case	4038917004770

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

This document was prepared by EDP support and is valid without signature.



SARSTEDT