

# Product specification

## Microvette® 200 K3 EDTA, 200 µl, cap violet, flat base



### Product information

Order number	20.1288.100
--------------	-------------

### Product descriptions

Name	Microvette® 200 K3 EDTA, 200 µl, Cap violet, Flat base
------	--

Product number	20.1288.100
----------------	-------------

Product description	Microvette® 200 K3 EDTA, capillary blood collection, preparation: K3 EDTA, nominal volume: 200 µl, (LxØ) with cap: 47.6 x 10.8 mm, with print, cap violet, flat base, inner receptacle: cylindrical, 100 piece(s)/case
---------------------	--

### Product characteristics

Base shape	Flat base
------------	-----------

Label/Print	With print
-------------	------------

Color of print/label	violet
----------------------	--------

Shape of inner tube	cylindrical
---------------------	-------------

Cap	Quadruple-threaded cap
-----	------------------------

Type of collection	capillary
--------------------	-----------

### Size

Diameter	10.8 mm
----------	---------

Length excluding cap	46.6 mm
----------------------	---------

Length including cap	47.6 mm
----------------------	---------

Sample volume	200 µl
---------------	--------

### Material & colours

Product material	Polypropylene (PP)
------------------	--------------------

Color of product	transparent
------------------	-------------

Cap material	High Density Polyethylene (HD-PE)
--------------	-----------------------------------

Color of cap

violet

## Purity & certification

Accordance with MDD or IVDD	In-vitro diagnostic, CE
CE certified	CE
Lot or serial number	yes

## Packaging

Piece(s) / inner box	100
Piece(s) / case	1000
Piece(s) / pallet	121000
Height of inner box	66 mm
Width of inner box	108 mm
Depth of inner box	139 mm
Height of case	160 mm
Width of case	228 mm
Depth of case	338 mm
Volume of case	0.0123 cbm
Weight of case	3.27 kg
EAN inner box	4038917127523
EAN case	4038917070393

**This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.**

**This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.**

This document was prepared by EDP support and is valid without signature.

