



Product specification

Microvette® 300 Fluoride/heparin FH, 300 µl, cap grey, flat base



Product information

Product number	20.1307.100
----------------	-------------

Product descriptions

Name	Microvette® 300 Fluoride/heparin FH, 300 µl, cap grey, flat base
------	--

Product number	20.1307.100
----------------	-------------

Product description	Microvette® 300 Fluoride/heparin FH, capillary blood collection, preparation: Fluoride + heparin, nominal volume: 300 µl, (LxØ) with cap: 47.6 x 10.8 mm, with print, cap grey, flat base, inner receptacle: conical, 100 piece(s)/case
---------------------	---

Material & colours

Color of product	transparent
Closure material	High Density Polyethylene (HD-PE)
Color of cap	grey

Size

Diameter	10.8 mm
Length excluding cap	46.6 mm
Length including cap	47.6 mm
Sample volume	300 µl

Product characteristics

Base shape	Flat base
Label/Print	With print
Color of print/label	grey
Shape of inner tube	conical
Cap	cap
Type of collection	capillary
Centrifugation recommended	2000 U/min
Centrifugation time	5 Min

Purity & certification

Product category	In-vitro diagnostic (IVD)
CE certified	CE - manufacturer's self-declaration
Batched	yes

Packaging

Piece(s) / inner box	100
Piece(s) / case	1000
Piece(s) / pallet	121000
Height of inner box	66 mm
Width of inner box	108 mm
Depth of inner box	139 mm
Height of case	160 mm
Width of case	228 mm
Depth of case	338 mm
Volume of case	0.0123 cbm
Weight of case	2.43 kg
EAN inner box	4038917127592
EAN case	4038917021272

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

