



Product specification

Microvette® 500 Serum CAT, 500 µl, cap red, flat base



Product information

Product number	20.1343.100
----------------	-------------

Product descriptions

Name	Microvette® 500 Serum CAT, 500 µl, cap red, flat base
Product number	20.1343.100
Product description	Microvette® 500 Serum CAT, capillary blood collection, preparation: clotting activator, nominal volume: 500 µl, (LxØ) with cap: 47.6 x 10.8 mm, with print, cap red, flat base, inner receptacle: cylindrical, 100 piece(s)/case

Material & colours

Product material	Polypropylene (PP)
Color of product	transparent
Closure material	High Density Polyethylene (HD-PE)
Color of cap	red

Size

Diameter	10.8 mm
Length excluding cap	46.6 mm
Length including cap	47.6 mm
Sample volume	500 µl

Product characteristics

Base shape	Flat base
Label/Print	With print
Color of print/label	red
Shape of inner tube	cylindrical
Cap	Quadruple-threaded cap
Type of collection	capillary
Centrifugation recommended	10000 x g
Centrifugation time	5 Min

Purity & certification

Product category	diagnostic device, CE
CE certified	CE
Batched	yes

Packaging

Minimum order quantity	1000
Piece(s) / inner box	100
Piece(s) / case	1000
Piece(s) / pallet	121000
Height of inner box	66 mm
Width of inner box	108 mm
Depth of inner box	139 mm
Height of case	160 mm
Width of case	228 mm
Depth of case	338 mm
Volume of case	0.0123 cbm
Weight of case	2.49 kg
EAN inner box	4038917127837
EAN case	4038917021425

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

This document was prepared by EDP support
and is valid without signature.



SARSTEDT