



## Product specification

### Microvette® 500 Serum Gel CAT, 500 µl, cap brown, flat base



## Product information

Product number	20.1344
----------------	---------

## Product descriptions

Name	Microvette® 500 Serum Gel CAT, 500 µl, cap brown, flat base
------	---

Product number	20.1344
----------------	---------

Product description	Microvette® 500 Serum Gel CAT, capillary blood collection, preparation: clotting activator/gel, nominal volume: 500 µl, (LxØ) with cap: 47.6 x 10.8 mm, with print, cap brown, flat base, inner receptacle: cylindrical, 100 piece(s)/case
---------------------	--

## Material & colours

Product material	Polypropylene (PP)
Color of product	transparent
Closure material	High Density Polyethylene (HD-PE)
Color of cap	brown

## Size

Diameter	10.8 mm
Length excluding cap	46.6 mm
Length including cap	47.6 mm

Sample volume	500 µl
---------------	--------

## Product characteristics

Base shape	Flat base
Label/Print	With print
Color of print/label	brown
Shape of inner tube	cylindrical
Cap	cap
Type of collection	capillary
Centrifugation recommended	10000 x g
Centrifugation time	5 Min

## Purity & certification

Product category	diagnostic device, CE
CE certified	CE
Batched	yes

## Packaging

Minimum order quantity	1000
Piece(s) / inner box	100
Piece(s) / case	1000
Piece(s) / pallet	121000
Height of inner box	66 mm
Width of inner box	108 mm
Depth of inner box	139 mm
Height of case	160 mm
Width of case	228 mm
Depth of case	338 mm
Volume of case	0.0123 cbm
Weight of case	2.66 kg
EAN inner box	4038917086264
EAN case	4038917004848

This document was prepared by EDP support  
and is valid without signature.



# SARSTEDT