



Product specification

Sample tube, Serum CAT, 4.5 ml, cap red, (LxØ): 75 x 13 mm, with paper label



Product information

Product number	32.329.100
----------------	------------

Product descriptions

Name	Sample tube, Serum CAT, 4.5 ml, cap red, (LxØ): 75 x 13 mm, with paper label
Product number	32.329.100
Product description	Sample tube, nominal volume: 4.5 ml, cap red, preparation: Clotting activator, push cap, (LxØ): 75 x 13 mm, round base, with paper label, 100 piece(s)/bag

Material & colours

Color of product	transparent
Tube material	Polypropylene (PP)
Closure material	Low Density Polyethylene (LD-PE)
Color of cap	red

Size

Diameter	13 mm
Length excluding cap	75 mm
Length of product	75 mm
Sample volume	4.5 ml

Product characteristics

Base shape	Round base
Label/Print	With paper label
Color of print/label	red
Graduation	Fill mark
Closure type	Push cap
Centrifugation recommended	2000 x g
Centrifugation time	10 Min

Purity & certification

Product category	In-vitro diagnostic (IVD)
CE certified	CE - manufacturer's self-declaration
Batched	yes

Packaging

Piece(s) / inner box	100
Piece(s) / case	1000
Piece(s) / pallet	64000
Height of case	215 mm
Width of case	298 mm
Depth of case	398 mm
Volume of case	0.0255 cbm
Weight of case	4.31 kg
EAN inner box	4038917105392
EAN case	4038917010542

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

This document was prepared by EDP support and is valid without signature.



SARSTEDT