



## Product specification

**Sample tube, Lithium heparin LH, 4 ml, cap green, (LxØ):  
75 x 12 mm, with paper label**



## Product information

Product number	36.336.100
----------------	------------

## Product descriptions

Name	Sample tube, Lithium heparin LH, 4 ml, cap green, (LxØ): 75 x 12 mm, with paper label
Product number	36.336.100
Product description	Sample tube, nominal volume: 4 ml, cap green, preparation: Lithium heparin, push cap, (LxØ): 75 x 12 mm, round base, with paper label, 100 piece(s)/bag

## Material & colours

Color of product	transparent
Tube material	Polypropylene (PP)
Closure material	Low Density Polyethylene (LD-PE)
Color of cap	green

## Size

Diameter	12 mm
Length excluding cap	75 mm
Length of product	75 mm
Sample volume	4 ml

## Product characteristics

Base shape	Round base
Label/Print	With paper label
Color of print/label	green
Graduation	Fill mark
Preparation concentration	16 IU/ml blood
Closure type	Push cap
Centrifugation recommended	2000 x g
Centrifugation time	10 Min

## Purity & certification

Product category	In-vitro diagnostic (IVD)
CE certified	CE - manufacturer's self-declaration
Batched	yes

## Packaging

Piece(s) / inner box	100
Piece(s) / case	2000
Piece(s) / pallet	60000
Height of case	340 mm
Width of case	398 mm
Depth of case	398 mm
Volume of case	0.0539 cbm
Weight of case	7.69 kg
EAN inner box	4038917105521
EAN case	4038917010641

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

This document was prepared by EDP support and is valid without signature.



# SARSTEDT