



## Product specification

**Sample tube, EDTA K3E, 4 ml, cap violet, (LxØ): 75 x 12 mm, with paper label**



## Product information

Product number	36.338.100
----------------	------------

## Product descriptions

Name	Sample tube, EDTA K3E, 4 ml, cap violet, (LxØ): 75 x 12 mm, with paper label
Product number	36.338.100
Product description	Sample tube, nominal volume: 4 ml, cap violet, preparation: K3 EDTA, push cap, (LxØ): 75 x 12 mm, round base, with paper label, 100 piece(s)/bag

## Material & colours

Color of product	transparent
Tube material	Polypropylene (PP)
Closure material	Low Density Polyethylene (LD-PE)
Color of cap	violet

## Size

Diameter	12 mm
Length excluding cap	75 mm
Length of product	75 mm
Sample volume	4 ml

## Product characteristics

Base shape	Round base
Label/Print	With paper label
Color of print/label	violet
Graduation	Fill mark
Preparation concentration	1.6 mg/ml blood
Cap	Skirted push cap
Closure type	Push cap

## Purity & certification

Product category	diagnostic device, CE
CE certified	CE
Batched	yes

## Packaging

Minimum order quantity	50000
Piece(s) / inner box	100
Piece(s) / case	2000
Piece(s) / pallet	60000
Height of case	340 mm
Width of case	398 mm
Depth of case	398 mm
Volume of case	0.0539 cbm
Weight of case	7.48 kg
EAN inner box	4038917105545
EAN case	4038917010658

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

This document was prepared by EDP support and is valid without signature.



# SARSTEDT