



Product specification

Micro sample tube Serum Gel CAT, 1.1 ml, screw cap, EU/ISO



Product information

Product number	41.1378.005
----------------	-------------

Product descriptions

Name	Micro sample tube Serum Gel CAT, 1.1 ml, screw cap, EU/ISO
Product number	41.1378.005
Product description	Micro sample tube Serum Gel CAT, venous blood collection, preparation: clotting activator/gel, nominal volume: 1.1 ml, (LxØ) with cap: 47 x 10.8 mm, screw cap, cap: brown, colour code EU/ISO, with paper label, label/print: brown, skirted conical base, 100 piece(s)/bag

Material & colours

Color of product	transparent
Tube material	Polypropylene (PP)
Closure material	Polypropylene (PP)
Color of cap	brown

Size

Diameter	10.8 mm
Length excluding cap	45 mm

Length including cap	47 mm
Length of product	47 mm
Sample volume	1.1 ml

Product characteristics

Base shape	skirted conical base
Label/Print	With paper label
Color of print/label	brown
Graduation	yes
Closure type	Screw cap
Type of collection	venous
Centrifugation recommended	10000 x g
Centrifugation time	5 Min

Purity & certification

Satisfies the requirement	IATA, ADR
Product category	In-vitro diagnostic (IVD)
CE certified	CE - manufacturer's self-declaration
Batched	yes

Packaging

Piece(s) / inner box	100
Piece(s) / case	1000
Piece(s) / pallet	96000
Height of case	140 mm
Width of case	278 mm
Depth of case	348 mm
Volume of case	0.0135 cbm
Weight of case	2.17 kg
EAN inner box	4038917127967
EAN case	4038917021500

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

This document was prepared by EDP support
and is valid without signature.



SARSTEDT