



## Product specification

### Micro sample tube Lithium heparin LH, 1.3 ml, screw cap, ISO



## Product information

Product number	41.1393.105
----------------	-------------

## Product descriptions

Name	Micro sample tube Lithium heparin LH, 1.3 ml, screw cap, ISO
------	--

Product number	41.1393.105
----------------	-------------

Product description	Micro sample tube Lithium heparin LH, venous blood collection, preparation: Lithium heparin, nominal volume: 1.3 ml, (LxØ) with cap: 47 x 10.8 mm, screw cap, cap: green, colour code ISO, with paper label, label/print: green, skirted conical base, 100 piece(s)/bag
---------------------	---

## Material & colours

Color of product	transparent
------------------	-------------

Tube material	Polypropylene (PP)
---------------	--------------------

Closure material	Polypropylene (PP)
------------------	--------------------

Color of cap	green
--------------	-------

## Size

Diameter	10.8 mm
----------	---------

Length excluding cap	45 mm
----------------------	-------

Length including cap	47 mm
Length of product	47 mm
Sample volume	1.3 ml

## Product characteristics

Base shape	skirted conical base
Label/Print	With paper label
Color of print/label	green
Graduation	yes
Preparation concentration	25 IU/ml blood
Closure type	Screw cap
Type of collection	venous
Centrifugation recommended	2000 x g
Centrifugation time	5 Min

## Purity & certification

Satisfies the requirement	IATA, ADR
Product category	In-vitro diagnostic (IVD)
CE certified	CE - manufacturer's self-declaration
Batched	yes

## Packaging

Piece(s) / inner box	100
Piece(s) / case	1000
Piece(s) / pallet	96000
Height of case	140 mm
Width of case	278 mm
Depth of case	348 mm
Volume of case	0.0135 cbm
Weight of case	1.93 kg
EAN inner box	4038917128322
EAN case	4038917021654

This document was prepared by EDP support  
and is valid without signature.



# SARSTEDT