



Product specification

Micro sample tube Lithium heparin LH, 1.3 ml, screw cap, ISO



Product information

Product number	41.1393.105
----------------	-------------

Product descriptions

Name	Micro sample tube Lithium heparin LH, 1.3 ml, screw cap, ISO
Product number	41.1393.105
Product description	Micro sample tube Lithium heparin LH, venous blood collection, preparation: lithium heparin LH, nominal volume: 1.3 ml, Ø: 10.8 mm, screw cap, cap: green, colour code ISO, with paper label, label/print: green, skirted conical base, 100 piece(s)/bag

Material & colours

Color of product	transparent
Tube material	Polypropylene (PP)
Closure material	Polypropylene (PP)
Color of cap	green

Size

Diameter	10.8 mm
Length excluding cap	45 mm
Length of product	47 mm

Sample volume	1.3 ml
---------------	--------

Product characteristics

Base shape	skirted conical base
Label/Print	With paper label
Color of print/label	green
Graduation	yes
Preparation concentration	25 IU/ml blood
Closure type	Screw cap
Type of collection	venous
Centrifugation recommended	2000 x g
Centrifugation time	5 Min

Purity & certification

Satisfies the requirement	IATA, ADR
Product category	diagnostic device, CE
CE certified	CE
Batched	yes

Packaging

Minimum order quantity	1000
Piece(s) / inner box	100
Piece(s) / case	1000
Piece(s) / pallet	96000
Height of case	140 mm
Width of case	278 mm
Depth of case	348 mm
Volume of case	0.0135 cbm
Weight of case	1.93 kg
EAN inner box	4038917128322
EAN case	4038917021654

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

This document was prepared by EDP support
and is valid without signature.



SARSTEDT