



## Product specification

### Micro sample tube Serum Gel CAT, 1.1 ml, push cap, EU/ISO



## Product information

Product number	41.1500.005
----------------	-------------

## Product descriptions

Name	Micro sample tube Serum Gel CAT, 1.1 ml, push cap, EU/ISO
Product number	41.1500.005
Product description	Micro sample tube Serum Gel CAT, venous blood collection, preparation: clotting activator/gel, nominal volume: 1.1 ml, (LxØ) with cap: 44 x 10.8 mm, push cap, cap: brown, colour code EU/ISO, with plastic label, label/print: brown, skirted conical base, 100 piece(s)/bag

## Material & colours

Color of product	transparent
Tube material	Polypropylene (PP)
Closure material	Low Density Polyethylene (LD-PE)
Color of cap	brown

## Size

Diameter	10.8 mm
Length excluding cap	40 mm

Length including cap	44 mm
Length of product	44 mm
Sample volume	1.1 ml

## Product characteristics

Base shape	skirted conical base
Label/Print	With plastic label
Color of print/label	brown
Graduation	yes
Cap	flat cap
Closure type	Push cap
Type of collection	venous
Centrifugation recommended	10000 U/min
Centrifugation time	5 Min

## Purity & certification

Product category	In-vitro diagnostic (IVD)
CE certified	CE - manufacturer's self-declaration
Batched	yes

## Packaging

Piece(s) / inner box	100
Piece(s) / case	1000
Piece(s) / pallet	96000
Height of case	140 mm
Width of case	278 mm
Depth of case	348 mm
Volume of case	0.0135 cbm
Weight of case	2.11 kg
EAN inner box	4038917129473
EAN case	4038917022248

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

This document was prepared by EDP support  
and is valid without signature.



# SARSTEDT