



Product specification

Micro sample tube Lithium heparin LH, 1.3 ml, push cap, ISO



Product information

Product number	41.1503.105
----------------	-------------

Product descriptions

Name	Micro sample tube Lithium heparin LH, 1.3 ml, push cap, ISO
Product number	41.1503.105
Product description	Micro sample tube Lithium heparin LH, venous blood collection, preparation: Lithium heparin LH, nominal volume: 1.3 ml, (LxØ) with cap: 44 x 10.8 mm, push cap, cap: green, colour code ISO, with plastic label, label/print: green, skirted conical base, 100 piece(s)/bag

Material & colours

Color of product	transparent
Tube material	Polypropylene (PP)
Closure material	Low Density Polyethylene (LD-PE)
Color of cap	green

Size

Diameter	10.8 mm
Length excluding cap	40 mm

Length including cap	44 mm
Length of product	44 mm
Sample volume	1.3 ml

Product characteristics

Base shape	skirted conical base
Label/Print	With plastic label
Color of print/label	green
Graduation	yes
Preparation concentration	25 IU/ml blood
Cap	flat cap
Closure type	Push cap
Type of collection	venous
Centrifugation recommended	2000 x g
Centrifugation time	5 Min

Purity & certification

Product category	diagnostic device, CE
CE certified	CE
Batched	yes

Packaging

Minimum order quantity	1000
Piece(s) / inner box	100
Piece(s) / case	1000
Piece(s) / pallet	96000
Height of case	140 mm
Width of case	278 mm
Depth of case	348 mm
Volume of case	0.0135 cbm
Weight of case	1.88 kg
EAN inner box	4038917129541
EAN case	4038917022286

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

This document was prepared by EDP support and is valid without signature.



SARSTEDT