



Product specification

Salivette®, with cotton swab, cap: white, with paper label



Product information

Product number	51.1534
----------------	---------

Product descriptions

Name	Salivette®, with cotton swab, cap: white, with paper label
Product number	51.1534
Product description	Salivette®, with cotton swab, cap: white, with paper label, label/print: white/black, for Saliva examination, 100 piece(s)/bag

Material & colours

Product material	Polypropylene (PP)
Closure material	Low Density Polyethylene (LD-PE)
Color of cap	white

Size

Diameter	16.8 mm
Length excluding cap	84 mm
Length including cap	97 mm

Product characteristics

Design	with cotton swab
Label/Print	With paper label
Color of print/label	white/black

Closure type	Push cap
Centrifugation recommended	1000 x g
Centrifugation time	2 Min

Purity & certification

Product category	diagnostic device, CE
CE certified	CE
Batched	yes

Packaging

Minimum order quantity	500
Piece(s) / inner box	100
Piece(s) / case	500
Piece(s) / pallet	28000
Height of case	125 mm
Width of case	398 mm
Depth of case	598 mm
Volume of case	0.0298 cbm
Weight of case	4.09 kg
EAN inner box	4038917086905
EAN case	4038917005432

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

This document was prepared by EDP support and is valid without signature.



SARSTEDT