



Product specification

Seraplas® valve filter, green, for separation of serum/plasma from the blood cells after centrifugat



Product information

Product number	53.420
----------------	--------

Product descriptions

Name	Seraplas® valve filter, green, for separation of serum/plasma from the blood cells after centrifugation
Product number	53.420
Product description	Seraplas® valve filter, V16, green, for separation of serum/plasma from the blood cells after centrifugation, 100 piece(s)/case

Material & colours

Color of product	green
------------------	-------

Size

Length of product	92 mm
-------------------	-------

Purity & certification

Product category	not regulated
CE certified	no CE mark
Batched	no
Safety data sheet	no

Packaging

Piece(s) / inner box	100
Piece(s) / case	1000
Piece(s) / pallet	54000
Height of inner box	130 mm
Width of inner box	99 mm
Depth of inner box	156 mm
Height of case	150 mm
Width of case	328 mm
Depth of case	513 mm
Volume of case	0.0252 cbm
Weight of case	2.98 kg
EAN inner box	4038917081498
EAN case	4038917000581

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

This document was prepared by EDP support and is valid without signature.



SARSTEDT